

Ultra Low Temperature Freezer



(LCD touch screen with USB port)

Product Feature

- Temperature Control

Microcomputer digital control. Temperature realization display. Display accuracy is 0.1℃

The temperature can be adjusted from -40 to -86 degrees.

Alarm Temperature can be set.

- Security System

With a perfect alarm system. Showing with sound and light .

Protection: Password protection, frequent start-up protection, etc.

Alarm types: High and low temperature alarm, sensor failure alarm

- Refrigeration System

The SECOP compressor make the performance stable.

The condensing fan brand in EBM,which made the noise in low.

Ultra-thick foam layer, double sealing of the door body, effective protection of the temperature inside the freezer.

100% CFC free.

- Ergonomic design

Easy Access: Built-in 304 stainless steel shelves can be adjusted.

High Security: The door is equipped with a security lock.

High Applicability: 187V-242V wide voltage design, suitable for unstable voltage environment.

- Certifications

CE, FDA,ISO9001,ISO13485

Specification

Model	UFZ-86V600
Capacity/L	608
Temperature	-40~-86℃
Power Supply/V/Hz	220/50
Refrigerant	mixed with CFC free
Noise /dB(A)	60
Compressor Brand	SECOP
Compressor Quantity	2
Shelves	3
Doors	4
Lock	YES
Wheels	YES
Alarm Way	Sound+Light
Alarm Type	High&low temperature alarm/Sensor failure alarm
Packing	
USB	Optional
Rated Input(w)	1100
Power Consumption	13.8 KWh/24H
Thermostat	Digital
Display	LED
Fan of Condenser	EBM
Cooling Way	Direct Cooling
Inner Size/mm	700*680*1255
Outer Size/mm	990*1000*1860
Packing Size/mm	1020*1080*2070
Foaming thickness/mm	120
Inner Material	304 stainless steel
Weight/kgs	265
Defrosting Way	Manual
Warranty	"1 year for machina 3 years for compressor"
LCD Screen	Optional