

Digital Ultrasonic Cleaner



Digital Ultrasonic Cleaner



UC-D



UC-D1/D2



UC-D3/D5



UC-D4/D6



UC-D7/D10/D15/D20

Digital Ultrasonic Cleaner

Advantage

- 1.High-power transducer for superior cleaning results
- 2.Digital LED display shows precise time and temperature.
- 3.Time setting 180s/360s or 1-30mins, built-in heating up to 80℃.
- 4.Super waterproof performance, longer service life
- 5.Main body, lid and base are made of high quality stainless steel
- 6.Use tap water, alcohol, organic solvents as cleaning fluid.



UC-D25



UC-D35/D40

Specification

Model	UC-D	UC-D1	UC-D2	UC-D3	UC-D5	UC-D4	UC-D6	UC-D7	UC-D10	UC-D15	UC-D20	UC-D25	UC-D30	UC-D45
Tank Capacity(L)	0.6	1.3	2	3.2	4.8	4.5	6.5	7	10	15	20	22	30	45
Transducer Qty(Pcs)	1	1	1	2	2	3	3	4	4	6	6	8	10	12
Ultrasonic Freq(KHz)	40													
Temperature Setting(℃)-Max	/ 80													
Drain Valve	No							Yes						
Time setting(min)	180S/360S	1-30 mins												
Cleaning Basket	Yes													
Carrying Handle	No			Yes										
Ultrasonic Power(W)-Max	50	70	70	120	120	180	180	240	240	360	360	480	600	720
Heating consumption(W)-Max	/	100	100	100	100	200	200	500	500	500	500	500	500	1000
Chamber Size (W×D×H)(mm)	150x85x65	150x135x65	150x135x100	240x135x100	240x135x150	300x150x100	300x150x150	300x240x100	300x240x150	330x300x150	330x300x200	500x300x150	500x300x200	500x300x300
Sizes (W×D×H)(mm)	175x110x155	180x165x205	180x165x240	265x165x235	265x165x285	330x180x260	330x180x310	330x270x260	330x270x310	360x330x330	360x330x380	530x330x340	530x330x390	530x330x490
Shipping size (W×D×H)(mm)	560x510x460	520x510x310	520x510x335	360x260x320	735x520x400	570x440x370	570x440x410	470×350×350	470×350×350	500x410×410	500×400×455	690×410×420	690×410×470	690×410×590
Packing Quantity	12pcs	4pcs	4pcs	1pc	4pcs	2pcs	2pcs	1pc	1pc	1pc	1pc	1pc	1pc	1pc
N.W./G.W. (kg)	1.6/22	2.6/13	2.8/14	3/4.5	4.1/21	5.2/14	5.8/25	7/9	7.8/10	10/13	10.5/13	13/16	15/18	19/23