

Bacti-cinerator Sterilizer

IS800 Series



Description:

- ◆ Bacti-cinerator Sterilizer is used in degassing for inoculating loops and needles, instead of alcohol lamps. Due to its convenience, sterilizer is widely used in the environment of biological safety cabinet, the exhaust fan, purification table and mobile van.

Application:

Widely used in biosafety cabinets, clean bench, near exhaust fan, mobile vehicles, even outdoor and other other severe environment.

Features:

- ◆ Safety, convenient. Totally take the place of alcohol lamps.
- ◆ Only 5-7 seconds are needed, because of its highest temperature of 825°C.
- ◆ Smart, light, easy to handle.
- ◆ Can be used in anaerobic chamber.
- ◆ Podzolic organics prevent pollutions effectively.
- ◆ Smart technique makes it possible to use for a long time and much safer.

Parameters:

Project	Parameters	
Product model	IS800-A	IS800-B
Max. temp. in central zone	825°C ± 25°C	825°C ± 25°C
Max. outer diameter of sterilization goods	φ 15mm	φ 35mm
Length of heating area	140mm	100mm
Heating time	15min	15min
Power	150W	250W
Voltage	AC110V or 220V, 50/60Hz	AC110V or 220V, 50/60Hz
Fuse	250V 3A 5X20	250V 3A 5X20
Ambient temperature	5°C ~ 40°C	5°C ~ 40°C
Dimensions(W×D×H)	120×155×180	130×120×200
Weight(kg)	1.0	1.3



Hangzhou Ruicheng Instrument Co., Ltd

ADD: Building E, No. 8 of Houshan Road, Yuhang District, Hangzhou City, Zhejiang Province, China
 TEL: 86-571-88630940 FAX: 86-571-85126370 Email: sales2@ruiwell.com
<http://www.ruiwell.net>