

Double Door/Three Door Refrigerator (2-8 ℃)



Double Door/Three Door Refrigerator (2-8 ℃)



Double Door/Three Door Refrigerator (2-8 ℃)

Specification

Model		MR-5V660	MR-5V1000	MR-5V1500
Capacity (L)		656	1006	1500
Size (WxDxH)	Internal (mm)	1100×454×1325	1100×684×1325	1680×595×1312
	External (mm)	1220×630×1885	1220×860×1885	1800×775×1965
Temperature	Range (℃) at RT.10~32 ℃	2~8(Optional 2~10; 2~14)	2~8(Optional 2~10; 2~14)	2~8(Optional 2~10; 2~14)
	Sensor	NTC	NTC	NTC
	Controller	Microprocessor	Microprocessor	Microprocessor
Cooling System	Refrigerant	R290 (CFC Free)	R290 (CFC Free)	R134a (CFC Free)
	Defrost	Auto	Auto	Auto
	Compressor Qty	1	1	Secop/1
Display		Digital display	Digital display	Digital display
Audible & Visible Alarm	Standard	High/low temperature, power failure, thermostat failure, door ajar	High/low temperature, power failure, thermostat failure, door ajar	High/low temperature, power failure, thermostat failure, door ajar
	Power Failure Backup System (Alarm)	8h	8h	8h
Material	Internal	Electro-galvanized steel with antimicrobial powder coating	Electro-galvanized steel with antimicrobial powder coating	304 stainless steel
	External	Electro-galvanized steel with antimicrobial powder coating	Electro-galvanized steel with antimicrobial powder coating	304 stainless steel
Accessories	USB Port	Optional	Optional	Optional
	Castors	Standard	Standard	Standard
	Chart Recorder	Optional	Optional	Optional
	Data Logger	Optional	Optional	Optional
	Remote Alarm Port	Optional	Optional	Optional
	Test Hole (Qty) / Φ (mm)	1/25	1/25	Optional
Shelves(Qty)		10 Shelves	10 Shelves	12 Shelves
Sound Emission (dB)		52	53	55
Electricity	Power Supply	220/50Hz/60Hz, 110/60Hz	220/50Hz/60Hz, 110/60Hz	220/50Hz/60Hz, 110/60Hz
	Consumption (W)	403	400	1065
N.W./G.W. (Kg)		166/184	189/213	245/280
Shipping Size (W×D×H) (mm)		1280×670×1950	1280×900×1950	1940×880×2035