

Mechanical Pipette

Eight-channel, Adjustable Volume,MP-EA series



Mechanical Pipette

Specification

MP-EA 8-channel Adjustable Volume Pipettes						
Volume Range	Increment	Test Volume	(Accuracy error)		(Precision error)	
			%	µL	%	µL
0.5-10µL	0.1µL	10µL	0.015	0.15	0.015	0.15
		5µL	0.025	0.125	0.025	0.125
		1µL	0.04	0.04	0.04	0.04
5-50µL	0.5µL	50µL	0.01	0.5	0.005	0.25
		25µL	0.015	0.375	0.01	0.25
		5µL	0.03	0.15	0.02	0.1
50-300µL	5µL	300µL	0.007	2.1	0.0025	0.75
		150µL	0.01	1.5	0.005	0.75
		50µL	0.015	0.75	0.008	0.4

Mechanical Pipette

Twelve-channel, Adjustable Volume,MP-TA series



Mechanical Pipette

Specification

MP-TA 12-channel Adjustable Volume Pipettes						
Volume Range	Increment	Test Volume	(Accuracy error)		(Precision error)	
			%	µL	%	µL
0.5-10µL	0.1µL	10µL	0.015	0.15	0.015	0.15
		5µL	0.025	0.125	0.025	0.125
		1µL	0.04	0.04	0.04	0.04
		50µL	0.01	0.5	0.005	0.25
5-50µL	0.5µL	25µL	0.015	0.375	0.01	0.25
		5µL	0.03	0.15	0.02	0.1
		300µL	0.007	2.1	0.0025	0.75
50-300µL	5µL	150µL	0.01	1.5	0.005	0.75
		50µL	0.015	0.75	0.008	0.4