

Multifunctional Ultrasonic Cleaner



Multifunctional Ultrasonic Cleaner



UC-M1



UC-M10



UC-M15

Advantage

- 1.Power and time are adjustable, temperature up to 60℃, time setting 1-99 mins.
- 2.Ultra-low noise
- 3.High-definition display shows accurate time, power and temperature.
- 4.Ergonomically designed carrying handle for easy mobility
- 5.Degassing and sweeping function
- 6.Long standby time with optional sleep mode for more energy saving.
- 7.Memory function can automatically save historical parameters.

Specification

Model	UC-M1	UC-M2	UC-M3	UC-M5	UC-M4	UC-M6
Tank Capacity(L)	1.3	2	3.2	4.8	4.5	6.5
Transducer Qty(Pcs)	1	1	2	2	3	3
Ultrasonic Freq(KHz)	40					
Temperature setting(℃)-Max	60℃					
Drain Valve	No	No	Yes			
Time setting(min)	1-99 mins					
Cleaning Basket	Yes					
Carrying Handle	No	No	Yes			
Power Supply	AC 220/110V,50/60Hz					
Ultrasonic Power(W)-Max	70(10%-100% adjustable)		120(10%-100% adjustable)		180(10%-100% adjustable)	
Heating consumption(W)-Max	100				200	
Chamber Size (W×D×H)(mm)	150x135x65	150x135x100	240x135x100	240x135x150	300x150x100	300x150x150
Sizes (W×D×H)(mm)	180x165x205	180x165x240	265x165x235	265x165x285	330x180x260	330x180x310
Shipping size (W×D×H)(mm)	520x510x310	520x510x335	730x520x340	730x520x400	430x350x290	420×260×390
Packing Quantity	4pcs	4pcs	4pcs	4pcs	1pc	1pc
N.W./G.W. (kg)	2.6/13	2.8/14	3.8/19	4.1/21	4.6/6	4.6/7



UC-M1



UC-M10/15

Multifunctional Ultrasonic Cleaner

Specification

Model	UC-M7	UC-M10	UC-M15	UC-M20	UC-M25	UC-M30	UC-M45
Tank Capacity(L)	7	10	15	20	22	30	45
Transducer Qty(Pcs)	4	4	6	6	8	10	12
Ultrasonic Freq(KHz)	40						
Temperature setting(℃)-Max	60℃						
Drain Valve	Yes						
Time setting(min)	1-99 mins						
Cleaning Basket	Yes						
Carrying Handle	Yes						
Power Supply	AC 220/110V,50/60Hz						
Ultrasonic Power(W)-Max	240(10%-100% adjustable)		360(10%-100% adjustable)		480(10%-100% adjustable)	600(10%-100% adjustable)	720(10%-100% adjustable)
Heating consumption(W)-Max	500						1000
Chamber Size (W×D×H)(mm)	300x240x100	300x240x150	330x300x150	330x300x200	500x300x150	500x300x200	500x300x300
Sizes (W×D×H)(mm)	330x270x260	330x270x310	360x330x330	360x330x380	530x330x340	530x330x390	530x330x490
Shipping size (W×D×H)(mm)	470x350x340	470x350x380	500×400×410	500×400×455	690×410×420	690×410×470	690×410×590
Packing Quantity	1pc	1pc	1pc	1pc	1pc	1pc	1pc
N.W./G.W. (kg)	7/9	7.8/10	10/13	10.5/13	13/16	15/18	19/23